

LabVIEW for Students Installation Instructions


Tuesday Mueller-Harder - 2024-06-10 - Comments (0) - Installation Instructions

Windows

1. Download the LabVIEW installer for Windows.
2. Right-click the installer and select **Run as Administrator**.
3. Select **Start** to unzip installation files.
4. Click on all components of LabVIEW you wish to install then click **Next**.
5. Click on any additional items you wish to install then click **Next**.
6. "Accept 2 Licenses", then click **Next**.
7. Click **Next**.
8. "Install Summary", **Next**.
9. Reboot to finish installer.


To Activate License in Windows: (after install)

1. Click Start Menu and launch "NI License Manager"
2. Click "Activate" Button on the top
3. Select "Enter serial number" in the dropdown



Use this license number for student use: **H12X04279**

***PLEASE NOTE: This above license works for all but 3 components of LabVIEW. There were be 3 remaining items that do not license with the Brown Student License.**



macOS

1. Download the LabVIEW disk image for macOS and double-click it.
2. From the LabVIEW volume, double-click the package file to install the base application.
3. Follow the installation steps, choosing the default options.
4. From the same volume, open the Extras folder and install the LabVIEW Runtime Engine using default choices.
5. Install any necessary drivers for your hardware from the Drivers page linked in the LabVIEW volume.
6. To install the Control Design and Simulation package, download and run the

LabVIEW_2022_Control_Design_and_Simulation_Module.pkg installer. Choose the default installation options.