

## Logger Pro Installation Instructions

Tuesday Mueller-Harder - 2023-10-25 - Comments (0) - Installation Instructions

### Windows

1. Download and run the Logger Pro installer for Windows.
2. You will receive a warning to disconnect your USB cable. Disconnect (if applicable) and click OK.
3. Follow the on-screen prompts, accepting all default installation options.
4. You may select "No thanks" if you do not wish to upgrade Go Pro to the latest version.
5. Once the installation is complete, **update to the current version** of Logger Pro immediately by clicking on the **Update** link on the Logger Pro catalog page.

### macOS

1. Download and run the Logger Pro installer for macOS.
2. Follow the on-screen prompts, accepting the default installation options.
3. Once installed, launch Logger Pro from your Applications folder.
4. **Update to the current version** of Logger Pro immediately by clicking on the **Update** link on the Logger Pro catalog page.